

WHITE LABS® ANALYTICAL SERVICES

BIG QC DAY

- TEST KIT -

**STOP
BEFORE
PROCEEDING**

- 1. Please read all instructions before starting.**
- 2. Register samples following instructions. Samples will not be processed if not registered.**
- 3. Please return samples at your earliest convenience. The results will be completed within two weeks of sample receipt.**

YOUR KIT INCLUDES THE FOLLOWING ITEMS FOR COLLECTION AND SHIPMENT OF TWO BEER SAMPLES.

- (2) 250mL sterile sample containers with unique barcodes**
- Return shipping label (except international customers)**
- Return shipping envelope**



PURE YEAST & FERMENTATION

1

REGISTER SAMPLES ONLINE

- Log in to your White Labs Sample Analysis Customer Portal at <https://whitelabs.qbench.net>.
- If you did not receive an email invitation to create a login or have forgotten your portal password, please email analyticalab@whitelabs.com or call (858) 693-3441 ext. 7387.
- Once you are signed in, click “Register Big QC Day Samples Here” at the top of the page.
- Select the Order ID.

Order Requests
White Labs Tasting Room / Order Requests / Filtered

Show 10 entries

ID	Request Status	Date Received
26405	PENDING CHANGES	07/17/2019 10:30 AM

Showing 1 to 1 of 1 entries

- Insert sample name in the description field for each Sample ID.

Samples

Sample ID	Sample Type	Description	Comments
39141	Beer		+ Add
39142	Beer		+ Add

[Save Samples](#)

- Click “Save Samples.”
- Click “Submit Sample Information” at the top right.

[Submit Sample Information](#)

2 OBTAIN CLEAN SAMPLES

All samples must be shipped in the supplied numbered containers.

Obtaining a beer sample from a beer bottle:

- Have barcoded sample container nearby with lid loosened.
- Use isopropanol spray to clean the outside top of beer bottle and bottle opener.
- Open the beer and carefully pour into the sample container. Fill to the 200mL mark.
- Close lid securely.

Obtaining a beer sample from tank:

- Have a barcoded sample container nearby with lid loosened.
- Clean the sample port with a cotton swab soaked with isopropanol.
- Use isopropanol spray to clean outside of the sample port. Flame faucet with portable flame if possible.
- Let approximately 8oz of beer flow through tap before collecting it into the sample container. Fill to the 200mL mark.
- Close lid securely.
- Repeat cleaning procedure once finished taking the sample.

3 SEND SAMPLES TO WHITE LABS

- Put containers in the return shipping envelope and seal. Put the return shipping label on shipping envelope.
- Call FedEx for a pick-up or drop off at your nearest FedEx location.

4 GET YOUR RESULTS

There are dates throughout the year that we run this test. Please reference our website for current dates at whitelabs.com/bigqcday.

Once testing is complete, you will receive an email with the reports. Your results will also be available to view via the White Labs Sample Analysis Customer Portal at whitelabs.qbench.net



[WHITELABS.QBENCH.NET](http://whitelabs.qbench.net)

[WHITELABS.COM](http://whitelabs.com) | [YEASTMAN.COM](http://yeastman.com)